



derma+flex® vs competition comparison chart



	<i>derma+flex</i> ®	Dermabond™	Indermil®	LiquiBand®	LiquiBand® Optima
Composition	60% Octyl, 40% Butyl	100% Octyl	100% Butyl	100% Butyl	10% Octyl, 90% Butyl
Flexibility	+++	+++	+	+	++
Strength	>30	10-30N	5-20N	5-20N	>30N
Viscosity	Non-drip gel	Low	Very Low	Very Low	Very Low
Application	Single layer	Apply in layers	Drip-on	Drip-on	Single layer
Pain on Application	Minimal	Significant	Minimal	Minimal	Minimal
Unit Size	0.7ml	0.5g	0.5g	0.5g	0.5g
Storage	Room Temp.	Room Temp.	Refrigeration	Refrigeration	Room Temp.

P/M.ME0008.1211.V5