

**derma+flex**<sup>®</sup>  
gel tissue adhesive

# High Viscosity Tissue Adhesive Like No Other!



## Derma+Flex Features & Benefits:

### 2-Octyl and N-Butyl Formulation

- Optimal Flexibility with 2-Octyl while maintaining superior strength of the N-Butyl
- Resistant to cracking or breaking

### Painless Application: No Sting or Burn

- Does not use an initiator or activator in its applicator tip, literally "Taking the Sting out of Medical Tissue Adhesives"

### Multiple Application Tips

- Broad surface applicator (self-puncturing cap and foam surface) allows spreading of the adhesive with uniformity and greater coverage
- Precision tip applicator (self-puncturing cap with an elongation) allows for detailed application

\*As long as adhesive remains intact

### Quick Drying

- Polymerises (dries) quickly within 60 seconds when applied correctly

### Waterproof Formula

- Prevents water from reaching the wound
- Just think how good your patient will feel after being able to shower within hours of their procedure instead of waiting for several days

### Multiple Uses on Single Patient

- Can be used on multiple wound sites on the same patient as its applicator tip will not harden while additional procedures are performed

### Microbial Barrier\*

- Protects the affected area by preventing bacteria from penetrating the wound
- Greatly reduces the chance of infection

### High Viscosity

- Syrup-like Viscosity at ~ 300cps
- Does not run or drip when applied correctly
- The perfect choice for vertical surface applications

